





Forsina Dicom Converter 3D



Quick Manual



About Forsina

"FORSINA" **the biggest medical library in the world**, aims to assist universities, students, healthcare professionals, and even patients in receive quality and accurate medical information, and it provides different tools to easily view, study, and understand human anatomy for male and female using 3D, Virtual Reality, and real CT scan models.

We have the biggest medically accurate content written developed by our anatomy experts and doctors, offering more than 1000 histology images, 50000 landmarks, 20000 traceable blood vessel and nerve paths, 500000 CT images and over 5 million medical words.

Forsina Products

Forsina has created software that mainly focus on the anatomy of the human body and histology in addition to radiology and Dicom convertor for the purpose of learning and harnessing the best tools to facilitate and empower the learning and researching process

Our Products



Forsina Anatomy 3D

Explore the female and male human body comprehensive structure in a stunning 3D environment covering every anatomical structure with full medical information for each one

Forsina Anatomy VR

Explore the female and male human body comprehensive structure in an immersive Virtual reality environment covering every anatomical structure with full medical information for each one



Forsina Dicom Converter 3D

Convert your DICOM files into a 3D module to be used in your desktop or mobile, this will allow you to explore your DICOM images in 3D environment

Forsina Dicom Converter VR

Convert your DICOM files into a Virtual Reality module to be used using your VR tool kit, this will allow you to explore your DICOM images in one of a kind experience you have never had before



Forsina Radiology 3D

Learn and understand how to read CT medical images, control the body density and frontal planes (axial, Sagittal, coronal) in a stunning 3D module created from over 120000 CT images for your studies

Forsina Radiology VR

Learn and understand how to read CT medical images, control the body density and frontal planes (axial, Sagittal, coronal) in an immersive VR module created from over 1 Million CT images for your studies

Quick Guide



Figure 1
start screen

Button	Function
Open file	Open ".vti" file in 3D environment
Convert Dicom	Convert Dicom files to ".vti" file
Convert CTA	Convert angiographic images into Dicom
Leave Forsina	Close the application














Figure 2
Dicom menu

Num.	Icon	Name	Function															
1		Density Presets	Display several preset densities for the user to choose from to view the model in the desired density.															
2		Volume Control Tool	<p>Reveal sliding bars that can be used to manually choose a density according to the user's needs. The shadow and color levels can also be manipulated using this tool.</p>															
3	<div><div>Slider Crop Cube Crop</div></div>	Cut Tool	<p>Cut the model in three different sectional planes (X, Y, Z). The cut tool consists of two separate cropping tools:</p> <ol style="list-style-type: none">Slider Crop: Reveal three sliders that you can move to cut the model in the desired plane.Cube Crop: Reveal a dotted cube around the model. Each dot on the cube can be held and moved to cut the model as desired. <p>- When selecting cube crop a small menu will appears</p> <div><table><tr><td>Pickable</td><td><input checked="" type="checkbox"/></td><td><u>Pickable</u>: if checked you can pick any dot to control it</td></tr><tr><td>Visibility</td><td><input type="checkbox"/></td><td><u>Visibility</u>: show/hide the cube outline</td></tr><tr><td>Face Handles</td><td><input checked="" type="checkbox"/></td><td><u>Face Handles</u>: show/hide dots in the middle of cube</td></tr><tr><td>Edge Handles</td><td><input checked="" type="checkbox"/></td><td><u>Edge Handles</u>: show/hide dots in the edges of cube</td></tr><tr><td>Corner Handles</td><td><input checked="" type="checkbox"/></td><td><u>Corner Handles</u>: show/hide dots in the corners of cube</td></tr></table></div>	Pickable	<input checked="" type="checkbox"/>	<u>Pickable</u> : if checked you can pick any dot to control it	Visibility	<input type="checkbox"/>	<u>Visibility</u> : show/hide the cube outline	Face Handles	<input checked="" type="checkbox"/>	<u>Face Handles</u> : show/hide dots in the middle of cube	Edge Handles	<input checked="" type="checkbox"/>	<u>Edge Handles</u> : show/hide dots in the edges of cube	Corner Handles	<input checked="" type="checkbox"/>	<u>Corner Handles</u> : show/hide dots in the corners of cube
Pickable	<input checked="" type="checkbox"/>	<u>Pickable</u> : if checked you can pick any dot to control it																
Visibility	<input type="checkbox"/>	<u>Visibility</u> : show/hide the cube outline																
Face Handles	<input checked="" type="checkbox"/>	<u>Face Handles</u> : show/hide dots in the middle of cube																
Edge Handles	<input checked="" type="checkbox"/>	<u>Edge Handles</u> : show/hide dots in the edges of cube																
Corner Handles	<input checked="" type="checkbox"/>	<u>Corner Handles</u> : show/hide dots in the corners of cube																

! Cropping the model can be done by simply moving the slider in the direction you want to crop.

Num.	Icon	Name	Function
4	 	CT Planes Color Window	Show 2D CT scan images of the model manipulate the color of the image or its background  
! Activating the Slider Crop tool while the CT plans tools is active the action on the crop tool will be reflected on the 2D CT image in all 3 axes			
5		Screenshot	Take a screenshot of the current view.
6		Save Current Density	Save the current density the user has set. This tool is intended to save time when you need to close the application and work on it later. 
7		Reset Camera	Reset the camera back to its original position.
8		Reset Volume	Reset the density, color, and position of the model to its original setting when the image was loaded.
9		New Volume	Open a new '.vti' file.



↑ Shift



To move the model, press and hold the shift and the left mouse buttons.

Ctrl



To roll the model, press and hold the Ctrl and the left mouse buttons.



To rotate the model, press and hold the left mouse while rotating